



DID SABZ



# ***Rotary Microtome***

**DS 8402**



شرکت دید سبز

# Rotary Microtome

**DS 8402**

یک سال گارانتی / ده سال پشتیبانی فنی

## Technical Data:

<b>Operating Temperature Range:</b>	+10°C to +35°C
<b>Section Thickness Range:</b>	1.0 - 60.0 μm
<b>Section Thickness Settings:</b>	
	1.0 μm increments from 1.0 – 10.0 μm
	2.0 μm increments from 10.0 – 20.0 μm
	5.0 μm increments from 20.0 – 60.0 μm
<b>Specimen Feed:</b>	approx. 30mm
<b>The Specimen Retraction:</b>	200 μm (can be turned off manually)
<b>Dimensions (W x D x H):</b>	413 x 618 x 305 mm
<b>Weight:</b>	approx. 34 kg
<b>Specimen Orientation:</b>	
<b>Horizontal</b>	8°
<b>Vertical</b>	8°
<b>Angle of Rotation:</b>	+90
<b>Trimming Stages:</b>	25 μm, 60 μm

## مشخصات و ویژگی ها:

- دارای پیشروی خیلی سریع با استفاده از هندول
- کنترل حد پیشروی سریع جلو و عقب در حالت دستی با استفاده از کلاج
- دارای فید برگشت اتوماتیک به عقب پس از اتمام عملیات برش در هر مرحله که امکان خاموش کردن آن وجود دارد.
- (Automatic Retraction)
- دارای لقمه گیر قابل تعویض با قابلیت نصب گیره ثابت جهت قالب های پلاستیکی یکبار مصرف
- کنترل حد مجاز پیشروی فیدهای میکرون توسط هندول بزرگ
- دارای هندول بسیار نرم و سبک
- امکان Trimming جهت دسترسی سریع به سطح بافت
- راحتی عملیات اپراتوری و کاملاً بی صدا

**DID SABZ Co.**

ADD.: Azadi Blvd., Opposite To Customs House  
Ostad Tanha Ave., Unit 6, Urmia-Iran  
Postal Code: 57176-66593  
Tel-Fax: (+98) 44 32771214-32781173



www.didsabz-co.com  
E-mail: info@didsabz-co.com

شرکت دید سبز (سهامی خاص)

آدرس: ارومیه، پستلوار آزادی، رویروی گسمرکه  
خجایان استاد تنها واحد ششبه دست راست  
شرکت دید سبز کد پستی: ۵۷۱۷۶-۶۶۵۹۳  
تلفن و موبایل: ۳۲۷۷۱۲۱۴-۳۲۷۸۱۱۷۳ (۰۴۴)