



DID SABZ



Dispenser Plus

DS4LM

Dispenser Plus

DS 4LM




شرکت دید سبز

یک سال گارانتی / ده سال پشتیبانی فنی

Technical Data:

System Control :	Microcontroller
Dimensions (W x D x H):	300 x 520 x 290 mm
Weight:	approx. 23.5 kg
Voltage:	220 v 50 Hz
Paraffin Capacity:	Up to 4 Liters

- Online temperature display with 0.1°C accuracy for tank and Hot plate
- Equipped with an electric valve which controls the flow volume
- Delay timer which could be set from 0 to 120 hours for holidays
- The alarm lights which conform to each function of the device
- The electric valve which could be set from 0 to 300 seconds and also could be turned off
- Two totally separate temperature controlling environments
- Controlling the paraffin tank using the thermal sensors from 30°C to 80°C
- Controlling the Hot Plate temperature from 30°C to 80°C using the thermal sensors
- Ease of reading the temperature via the flat panel holding 7 segment display



In Conformity to IEC61010-1

مشخصات و ویژگی ها:

- سیستم کنترل دستگاه از نوع میکروکنترلر
- راحتی قرائت دمای پانل فلت با نمایش 7segment
- مجهز به پدال پایی برای آزاد بودن هر دو دست اپراتور جهت بالا بردن دقت و سرعت عمل قالب گیری
- دو محیط کنترل دما بصورت مجزا
- کنترل قطع و وصل سیال توسط شیر الکتریکی بدون چکه و با قطعی کامل
- جلوگیری از انسداد مسیر خروجی پارافین با تامین حرارت لازم
- کنترل مقدار حجم خروجی سیال با استفاده از یک سیستم الکترونیکی
- ظرفیت مخزن ۴ لیتر
- تایمر تاخیری قابل برنامه ریزی برای ایام تعطیلات از صفر تا ۱۲۰ ساعت معادل ۵ روز
- نصب چراغ های هشدار دهنده متناسب عملکرد دستگاه
- شیر برقی قابل کنترل توسط تایمر صفر تا ۳۰۰ ثانیه که می توان تایمر را خاموش کرد.
- کنترل صفحه گرم در محدوده دمایی از ۳۰ الی ۸۰ درجه سانتیگراد با استفاده از سنسور حرارتی

DID SABZ Co.

ADD: Azadi Blvd., Opposite To Customs House
Ostad Tanha Ave., Unit 6, Urmia-Iran
Postal Code: 57176-66593
Tel-Fax: (+98) 44 32771214-32781173



www.didsabz-co.com
E-mail: info@didsabz-co.com

شرکت دید سبز (سهامی خاص)

آدرس: ارومیه - پلوار آزادی روی سروی گسمرکه
خیابان استاد تنها واحد ششده دست راست
شرکت دید سبز کد پستی: ۵۷۱۷۶-۶۶۵۹۳
تلفن و موبایل: ۳۲۷۸۱۱۷۳-۳۲۷۷۱۲۱۴ (۰۴۴)